




## **Technical Data**



## **Tracoe Comfort Plus**

Size	ID- OC	OD-OC bc	OD-OC bf	ID- IC	TL	ОВ	ВА
mm	mm	mm	mm	mm	mm	mm	Θ°
07	7.0	8.8	10.5	5.8	81	87	100
08	8.0	10.2	11.3	6.8	84	91	100
09	9.0	11.6	12.6	7.6	86	94	100
10	10.0	12.0	13.7	8.7	88	97	100
11	11.0	13.5	15.0	9.7	91	100	100

ID-OC: inner diameter (clear width) at bottom of outer cannula;
OD-OC bc: outer diameter at bottom of outer cannula; OD-OC bf: outer
diameter of outer cannula behind the flange; ID-IC: inner diameter at bottom
of inner cannula; TL: length along center line from start of neck flange to
bottom of tube; OB: length along outer bend (greatest length); BA: bending
angle

Please be aware when switching a patient from another tube type in a certain size to one of the Tracoe Comfort tubes, the outer and inner diameter should be considered.

## Head office:

Atos Medical AB Hyllie Boulevard 17 SE-215 32 Malmo, Sweden Tel: +46 (0) 415 198 00 Email: info@atosmedical.com

## Manufacturer:

Tracoe Medical GmbH Reichelsheimer Straße 1/3 55268 Nieder-Olm, Germany MC3049-TcEN 202306

© Coloplast A/S. All rights reserved. Atos and the Atos Medical logo are trademarks of Coloplast A/S.