



Technical Data

Tracoe Vario

Size	ID	OD	OB	CD	BA	REF 450 / REF 450-P	REF 455	REF 460 / REF 460-P	REF 464	REF 470 / REF 470-P
						TL (min-max) mm	TL (min-max) mm	TL (min-max) mm	TL (min-max) mm	TL (min-max) mm
	mm	mm	mm	mm	∅°					
06	6.0	8.2	76	18.0	100	50.0 - 62.0	59.5 - 71.5	50.0 - 62.0		
07	7.0	9.7	86	23.0	100	53.0 - 68.0	64.5 - 79.5	54.0 - 69.0	64.0 - 79.0	55.0 - 70.0
08	8.0	11.2	97	28.0	100	64.0 - 84.0	71.5 - 91.5	65.0 - 85.0	71.0 - 91.0	63.0 - 83.0
09	9.0	12.3	106	30.0	100	66.0 - 93.0	76.0 - 103.0	64.0 - 91.0	76.0 - 103.0	64.0 - 91.0
10	10.0	13.7	120	32.0	100	63.0 - 98.0	75.0 - 110.0	68.5 - 103.5		64.0 - 99.0
11	11.0	14.8	121	34.0	100	68.0 - 103.0	77.0 - 112.0	68.0 - 103.0		

ID: inside diameter (clear width) at bottom of cannula; **OD:** outside diameter at bottom of cannula;
TL min: minimal length along centre line; **TL max:** maximal length along centre line; **OB:** length along outer bend;
BA: bending angle; **CD:** cuff diameter

Head office:
 Atos Medical AB
 Hyllie Boulevard 17
 SE-215 32 Malmö, Sweden
 Tel: +46 (0) 415 198 00
 Email: info@atosmedical.com

Manufacturer:
 Tracoe Medical GmbH
 Reichelsheimer Straße 1/3
 55268 Nieder-Olm, Germany

MC3045-TcEN 202306

© Coloplast A/S. All rights reserved.
 Atos and the Atos Medical logo are
 trademarks of Coloplast A/S.